



Dear Valued Customer,

In our continuing effort to address your practice's needs, we have implemented the following change. Starting July 1st, we will no longer stamp progressive digital lenses. Due to differences in base curves, parallax, etc., these stamps are often not properly applied, and are difficult to remove. All lenses will be marked right (R) and left (L) and all lenses will be dotted on the nasal and temple circles.


As always, we appreciate your business and look forward to hearing from you.


Progressive lenses

Marking changes

Before

1. Right (R) or Left (L) on top of the lens
2. Lens with nasal and temple dots
3. Yellow stamp






After

1. Right (R) or Left (L) on top of the lens
2. Lens with nasal and temple dots

⚠ No more progressive yellow stamp



Benefits of the change:

- ✓ Progressive stamp is not very accurate. It's better to use the dots as reference;
- ✓ It's not necessary to wash the lenses to remove the yellow ink;
- ✓ Reduce the risk of scratches.

